

Although the polishing process is imperfect, leaving scratches on the section surface, adequate optical properties are preserved. Most of the abrasions disappear when the liquid cover glass is applied, and the microscope is usually focused below those that remain. The time saved by imperfect polishing is substantial, and the risk of sample damage from prolonged handling is reduced by this rapid procedure.

References

1. Dahl, L. K.: A simple and sensitive histochemical method for calcium. *Proc. Soc. exp. Biol. (N.Y.)* **80**, 474-479 (1952)
2. Frost, H. M.: Tetracycline-based histological analysis of bone remodeling. *Calcif. Tiss. Res.* **3**, 211-237 (1969)
3. Galante, J., Rostoker, W., Lueck, R., Day, R. D.: Sintered fiber metal composites as a basis for attachment of implants to bone. *J. Bone Jt Surg.* **A53**, 101-114 (1971)
4. Huber, J. D., Parker, F., Odland, G. F.: A basic fuchsin and alkalinized methylene blue rapid stain for epoxy-embedded tissue. *Stain Technol.* **43**, 83-87 (1968)
5. Hulbert, S. F., Young, F. A., Mathews, R. S., Klawitter, J. J., Talbert, C. D., Stelling, F. H.: Potential of ceramic materials as permanently implantable skeletal prostheses. *J. Biomed. Mater. Res.* **4**, 433-456 (1970)
6. Karagianes, M. T.: Porous metals as a hard tissue substitute. *Biomat. Med. Devel. Artif. Organs* **1**, 171-181 (1973)
7. Morrison, S. J.: Tissue reaction to three ceramics of porous and non-porous structures. Masters Thesis, Div. Interdisciplinary Studies, Clemson University, Clemson, S.C., May 1971
8. Spurr, A. R.: A low-viscosity epoxy resin embedding medium for electron microscopy. *J. Ultrastruct. Res.* **26**, 31-43 (1969)

Erratum

E. Hausmann, R. Genco, N. Weinfeld and R. Sacco: "Effects of Sera on Bone Resorption in Tissue Culture". Vol. 13, pp. 311-317 (1973).

In Table 3, experiment number 84, A) rat treatment should have read "Ag-Ab ppt".