

MEETING ABSTRACT

Open Access

Quantification and accuracy of clinical [11C]-PiB PET/MRI: the effect of MR-based attenuation correction

Ian Law^{1*}, Flemming L Andersen¹, Adam E Hansen¹, Steen G Hasselbalch², Claes Ladefoged¹, Sune H Keller¹, Søren Holm¹, Liselotte Højgaard¹

From PSMR14: 3rd Conference in PET/MR and SPECT/MR
Kos Island, Greece. 19-21 May 2014

¹The Department of Clinical Physiology, Nuclear Medicine and PET, Rigshospitalet, Denmark

The Dixon-Water-Fat segmentation (DWFS) method is a standard attenuation correction (AC) method in PET/MRI on the Siemens mMR and has demonstrated a systematic quantitative bias in [18F]-FDG-PET/MRI studies of the brain compared to PET/CT. The aim of this study was to evaluate the impact of DWFS-AC in a hybrid PET/MR scanner on regional and global quantitation of [11C]-PiB cerebral amyloid imaging of the brain and in the clinical reading.

Twenty-eight healthy volunteers, 6 with mild cognitive impairment (MCI), 4 with Alzheimers disease (AD), and 6 other dementia cases underwent a simultaneous PET/MRI (Siemens mMR) acquisition 40 min pi of (170-709) MBq [11C]-PiB. A single 30 min frame was reconstructed twice (OSEM-3D, 2lit 4 sub, 5 mm Gauss) with the only difference being the AC calculation. AC was performed using either the standard DWFS or a head low-dose CT scan acquired independently on the same day. Activity concentration (Bq/mL) was sampled on AC-PET from symmetrically delineated ROI's. Ratios of region-to-cerebellar grey matter (SUVr) were calculated. Amyloid uptake was considered abnormal at cortical SUVr >1.5 or increased amyloid uptake in two or more grey matter regions on visual evaluation.

The average activity concentration in all ROI's was biased by -18 % in DWFS-AC compared to CT-AC. Although the visual categorization of both, amyloid-positive (n=11) and -negative (n=33) scans was unaffected by DWFS-AC, 3 healthy subjects were quantitatively reclassified as amyloid positive. The SUVr values were overestimated by 0.11 in the caudate nuclei and 0.04 in lateral cortical ROI's in DWFS-AC compared to CT-AC.

The visual and quantitative consequences of MR-AC using DWFS in the brain with [11C]-PiB exhibit a noticeable radially variable bias. Although robust in visual evaluation, the quantitative diagnostic criteria (SUVr) using this biomarker with DWFS-AC may need to be modified.

Authors' details

¹The Department of Clinical Physiology, Nuclear Medicine and PET, Rigshospitalet, Denmark. ²Memory Disorders Research Group, Rigshospitalet, Denmark.

Published: 29 July 2014

doi:10.1186/2197-7364-1-S1-A69

Cite this article as: Law et al.: Quantification and accuracy of clinical [11C]-PiB PET/MRI: the effect of MR-based attenuation correction. *EJNMMI Physics* 2014 **1**(Suppl 1):A69.

Submit your manuscript to a SpringerOpen[®] journal and benefit from:

- ▶ Convenient online submission
- ▶ Rigorous peer review
- ▶ Immediate publication on acceptance
- ▶ Open access: articles freely available online
- ▶ High visibility within the field
- ▶ Retaining the copyright to your article

Submit your next manuscript at ▶ springeropen.com
