

The Eyelink Toolbox: Eye tracking with MATLAB and the Psychophysics Toolbox

FRANS W. CORNELISSEN and ENNO M. PETERS
University of Groningen, Groningen, The Netherlands

and

JOHN PALMER
University of Washington, Seattle, Washington

The Eyelink Toolbox software supports the measurement of eye movements. The toolbox provides an interface between a high-level interpreted language (MATLAB), a visual display programming toolbox (Psychophysics Toolbox), and a video-based eyetracker (Eyelink). The Eyelink Toolbox enables experimenters to measure eye movements while simultaneously executing the stimulus presentation routines provided by the Psychophysics Toolbox. Example programs are included with the toolbox distribution. Information on the Eyelink Toolbox can be found at <http://psychtoolbox.org/>.

Measuring eye movements during psychophysical tasks and experiments is important for studying eye-movement control, gaining information about a level of behavior generally inaccessible to conscious introspection, examining information processing strategies, as well as controlling task performance during experiments that demand fixation or otherwise require precise knowledge of a subject's gaze direction (e.g., Brenner & Cornelissen, 2000; Cornelissen & Dobbelsteen, 1999). Eye-movement recording is becoming a standard part of psychophysical experimentation.

Although eye-tracking techniques exist that rely on measuring electrical potentials generated by the moving eye (electro-oculography) or a metal coil in a magnetic field, such methods are relatively cumbersome and uncomfortable for the subject (e.g., because electrodes have to be attached to the head or a coil has to be placed on the cornea). A new generation of eye trackers is now available based on the noninvasive recording of images of a subject's eye using infrared sensitive video technology and relying on software processing to determine the position of the subject's eyes relative to the head. Since these trackers are video based, there is no need for direct contact with the subject's eyes, making these trackers much more suitable for routine eye-movement recording during longer sessions. By combin-

ing the information on eye position with a measure of head position, estimates of gaze position on a display can be obtained, allowing the creation of gaze-dependent displays.

Eyelink Gazetracker

The Eyelink Gazetracker (SR Research Ltd., Mississauga, Ontario, Canada) is one of these video-based eye trackers and is used in such research fields as psychology, ophthalmology, neurology, and ergonomics. The Eyelink uses two high-speed cameras (CCD sensors) to track both eyes simultaneously. A third camera (also a CCD sensor) tracks four infrared markers mounted on the visual stimulus display, so that head motion can be measured and gaze position can be computed. The cameras produce images at a sampling rate of 250 Hz (4-msec temporal resolution). The Eyelink is used with a PC with dedicated hardware for doing the image processing necessary to determine gaze position. Pupil position is tracked by an algorithm similar to a centroid calculation, with a noise-limited resolution of 0.01° or less, and a velocity noise of less than 3°/sec (manufacturer's specifications). An optional heuristic filter (Stampe, 1993) removes single-sample noise artifacts and does not affect measured saccade velocity and acceleration. Pupil position is mapped to a head-referenced coordinate system by a biquadratic mapping function (Sheena & Borah, 1981; Stampe, 1993) for two-dimensional eye tracking, or by a quadratic function for one-dimensional eye tracking. Head-position data are then combined with the head-referenced data to compute the true gaze position on the stimulus display. The optical head-tracking system has extremely low angular noise and therefore does not appreciably increase the final noise level of the system.

The Eyelink communicates via a high-speed ethernet connection with a separate computer (Macintosh or PC) that performs the stimulus display. The Eyelink system

Eyal Reingold and Dave Stampe conceived and developed the Eyelink Gazetracker. We greatly appreciate the effort David Brainard and Denis Pelli have undertaken in developing the Psychophysics Toolbox (see <http://psychtoolbox.org/>). We thank Eyal Reingold, Francesco Maringelli, Erin Harley, and two anonymous referees of an earlier version of this paper for commenting on the manuscript. J.P. was partly supported by HHMI and NIH Grant EY11378. Correspondence should be addressed to F. W. Cornelissen, Laboratory of Experimental Ophthalmology, School for Behavioral and Cognitive Neurosciences, University of Groningen, P.O. Box 30001, 9700 RB Groningen, The Netherlands (e-mail: f.w.cornelissen@med.rug.nl).

Table 1
Listing and Description of the Main Commands Available in the Eyelink Toolbox

Command	Function
'Initialize'	Establishes an ethernet connection between the Eyelink PC and display computer and initializes the Eyelink Gazetracker.
'Initializedummy'	Run in "dummy mode," which allows one to run (and, e.g., debug) stimulation programs without a Gazetracker being connected.
'Shutdown'	Close the connection to the Eyelink PC and Gazetracker.
'Isconnected'	Reports the status of the connection (connected, no connection, running in dummy mode).
'Openfile'	Opens a file on the Eyelink PC to store data and events in.
'Closefile'	Closes the data file.
'Command'	Sends a command string to the Eyelink PC (e.g., 'link_sample_data=GAZE' tells the Eyelink to send real-time gaze data to the display computer).
'Message'	Sends a message to the Eyelink PC (e.g., "Stimulus On").
'Initwindow'	Provide information about the graphics environment (e.g., screen resolution) to the Eyelink Toolbox.
'Eyeavailable'	Reports which eye(s) is being tracked.
'Trackersetup'	Perform integrated tracker set-up procedure.
'Dodriftcorrect'	Perform integrated drift correction procedure.
'Startrecording'	Starts recording eye-movement data.
'Stoprecording'	Stops recording eye-movement data.
'Checkrecording'	Reports whether the Eyelink is recording data.
'Newsampleavailable'	Report whether a new data sample (either in int or float format) is available (for real-time purposes).
'Newfloatsampleavailable'	
'Newestsample'	Return the most recent data sample (either in int or float format).
'Newestfloatsample'	

parses the eye-movement data on line and in real time so that information on eye position is available almost instantaneously. Average delays from eye movement to position data availability are 6 msec with heuristic filtering disabled, and 10 msec with filtering enabled (manufacturer's specifications). This allows the creation of (near) real-time gaze-contingent applications (e.g., Cornelissen & Dobbelsteen 1999; Tant, Cornelissen, Kooijman, & Brouwer, 2002). The parser also provides information on fixation and saccade events and associated statistics immediately after the events have been completed.

One limitation of the Eyelink is that there is considerable trial-to-trial variability in the estimate of absolute gaze position. For example, typical measurements of the standard deviations of gaze position at fixation are 0.5° horizontal and 1.5° vertical. This variability might be due to slippage of the headmounted system over time. This problem can be minimized by using fixation position at the beginning of each trial to "correct" the gaze estimate (*drift correction*). This method provides good relative position estimates within a trial at the cost of information about absolute position across trials. Saccade parameter and fixation position estimates of the Eyelink were highly correlated with those measured using a scleral coil system. The Eyelink's relatively low sampling frequency of 250 Hz results in somewhat noisier parameter estimates in case of small saccades (van der Geest & Frens, 2002).

Stimulus Presentation

Most eye-movement experiments require the display of visual targets for the subject to track or look at. Computer displays have virtually become standard equipment for conducting visual psychophysics experiments since these allow precise software specification of the stimu-

lus. Although low-level programming languages such as C enable the required accuracy in control of timing, luminance, and color output of displays, they do not provide a friendly programming environment. The Psychophysics Toolbox (Brainard, 1997; Pelli, 1997) is a software package that adds the ability for precise stimulus specification to MATLAB, a high-level interpreted language with extensive support for numerical calculations (The MathWorks, 1993), allowing for rapid and flexible programming of psychophysical experiments.

We have developed the Eyelink Toolbox, a high-level interface between MATLAB and the Eyelink Gazetracker. The toolbox enables one to measure eye movements while simultaneously executing stimulus presentation routines provided by the Psychophysics Toolbox as well as other MATLAB scripts. The powerful combination of the Psychophysics and Eyelink Toolboxes allows for a

Table 2
Listing of the Main Fields Available in the Eyelink Data Sample Structure

Field	Content
time	Time stamp of sample.
flags	Flags to define what data are included in each sample.
gx, gy	Horizontal and vertical gaze position data (x and y) in screen coordinates (pixels) for the left and right eyes (depending on whether they are actually being tracked).
pa	Pupil size or area (depending on settings of Eyelink) of left and right eyes.
rx, ry	Horizontal and vertical resolution of the gaze data (pixels per degree) for the left and right eyes.
hx, hy	Horizontal and vertical head-referenced eye-position data for the left and right eyes.

Note—Whether these fields are actually filled with useful information depends on the settings of the Eyelink. This can be changed at run time using the 'command' command of the Eyelink Toolbox (see Table 1).

relatively fast and easy implementation of experiments involving eye tracking (e.g., gaze-dependent displays).

Implementation of the Toolbox

The toolbox consists of a combination of a MATLAB extension (MEX) file and code written in native MATLAB language. The extension, which is partly based on proprietary low-level C subroutines provided by the manufacturer of the eye tracker, can be called directly from MATLAB. For timing and graphics manipulation, the Eyelink Toolbox relies on routines of the Psychophysics Toolbox (Brainard, 1997) and VideoToolbox (Pelli, 1997).

The Eyelink Toolbox provides access to all Eyelink routines. A listing of the main commands and a short description of their function is provided in Table 1. In addition, the Eyelink Toolbox provides both integrated routines for performing calibration and drift correction, as well as procedures written in native MATLAB language that allow for a large degree of customization (e.g., of calibration targets and routines; see the listing in the Appendix). To have (near) real-time gaze-dependent displays, the Eyelink parses the eye-movement data on line and sends gaze data, as well as other data, over the ethernet connection to the display computer. The Eyelink Toolbox allows access to this information via a MATLAB structure. Table 2 lists the main fields of this structure and provides a short description of its contents.

The use of an interpreted language comes at a cost: Native MATLAB code tends to run significantly slower than comparable C code. The penalty of the overhead is about 4 msec per MATLAB statement (as reported by the PsychToolbox's 'SpeedTest.m' routine). Nevertheless, the Eyelink Toolbox implementation, even on a relatively old Power Macintosh with a 400-MHz G3 processor, is fast enough to enable a near real-time gaze-dependent display. Between eye movement and screen update there is a delay of about 20 msec,¹ of which our tests indicate that almost all is due to the Eyelink hardware and software on the dedicated Eyelink PC. A millisecond or less of the delay is due to network communication between the Eyelink and display computers. Only a fraction of a millisecond is due to the overhead of using MATLAB in addition to the underlying C routines. Thus, the use of the Eyelink Toolbox within the MATLAB environment, instead of directly programming in C, adds little cost in terms of execution time. The MATLAB environment, by providing extensive support for numerical processing, in principle also allows the postprocessing of eye-movement data within the same environment. At present, support for this additional use of MATLAB is not part of the Eyelink Toolbox.

Example and Use

The Appendix lists an example of a short MATLAB program that uses the Eyelink and Psychophysics Toolboxes to create a simple gaze-dependent display. Step 1 in the example is the initialization of the Eyelink; Step 2 opens a graphics window using the Psychophysics Tool-

box SCREEN routine. In Step 3, the Eyelink Toolbox default values are set, and information about the graphics environment is passed onto the Eyelink mex routine, and, in Step 4, the Eyelink's integrated calibration and drift correction procedures are performed. In Step 5, data recording commences. In Step 6, current gaze position as communicated by the Eyelink is used to show a small dot on the display that moves with the subject's gaze. In case the Eyelink provides no valid gaze data (e.g., during a blink), the screen is immediately blanked. In Step 7, graphics window, data file, and tracker are closed.

Documentation and Availability

The Eyelink Toolbox can be downloaded at <http://psych-toolbox.org/> and, like the Psychophysics Toolbox, is available for both the Macintosh and the Windows operating systems.² Installation consists of adding a single folder (approximately 1 MB) to MATLAB's collection of toolboxes. Help is provided via MATLAB's regular conventions. The distribution includes example MATLAB code (e.g., the program listed in the Appendix). The C source code of the toolbox is available (though not for the proprietary subroutines on which the toolbox is based). Note that the Eyelink Toolbox is neither provided nor endorsed nor supported by the manufacturer of the Eyelink Gaze-tracker. The toolbox may be used freely for teaching or research. It may not be used for commercial gain without permission of the first author of the present article.

Conclusion

The Eyelink Toolbox, in combination with the Psychophysics Toolbox, provides a fast, easy, interactive, and powerful means to develop research-grade eye-movement paradigms (see, Huk, Palmer, & Shadlen, 2002; Li, Brenner, Cornelissen, & Kim, 2002).

REFERENCES

- BRAINARD, D. H. (1997). The Psychophysics Toolbox. *Spatial Vision*, **10**, 437-442.
- BRENNER, E., & CORNELISSEN, F. W. (2000). Separate simultaneous processing of egocentric and relative positions. *Vision Research*, **40**, 2557-2563.
- CORNELISSEN, F. W., & DOBBELSTEEN, J. (1999). Heading detection with simulated visual field defects. *Visual Impairment Research*, **1**, 71-84.
- CORNELISSEN, F. W., & KOOIJMAN, A. C. (2002). *The influence of artificial scotomas on eye-movements during visual search*. Manuscript submitted for publication.
- HUK, A. C., PALMER, J., & SHADLEN, M. N. (2002). Temporal integration of visual motion information: Evidence from response times. *Journal of Vision*, **2** (7), 228a.
- LI, H.-C. O., BRENNER, E., CORNELISSEN, F. W., & KIM, E. S. (2002). Systematic distortion of perceived 2D shape during smooth pursuit eye movements. *Vision Research*, **42**, 2569-2575.
- PELLI, D. G. (1997). The VideoToolbox software for visual psychophysics. *Spatial Vision*, **10**, 437-442.
- SHEENA, D., & BORAH, B. (1981). Compensation for some second-order effects to improve eye position measurements. In D. F. Fisher, R. A. Monty, & J. W. Senders (Eds.), *Eye movements: Cognition and visual perception*. Hillsdale, NJ: Erlbaum.
- STAMPE, D. M. (1993). Heuristic filtering and reliable calibration methods for video-based pupil-tracking systems. *Behavior Research Methods, Instruments, & Computers*, **25**, 137-142.

- TANT, M. L. M., CORNELISSEN, F. W., KOOIJMAN, A. C., & BROUWER, W. H. (2002). Hemianopic visual field defects elicit hemianopic scanning. *Vision Research*, **42**, 1339-1348.
- THE MATHWORKS. (1993). MATLAB user's guide. The MathWorks, Inc., Natick, MA.
- VAN DER GEEST, J. N., & FRENS, M. A. (2002). Recording eye movements with video-oculography and scleral search coils: A direct comparison of two methods. *Journal of Neuroscience Methods*, **114**, 185-195.

NOTES

1. We determined the delay from eye movement to stimulus update to be about 20 msec. This was determined by using electro-oculography to independently record eye movements and a photodiode to measure stimulus-update-related changes in light flux on the screen (Cornelissen & Kooijman, 2002).
2. The EYELINK Toolbox is compatible with both the first and the second generations of EYELINK Gazetrackers.

APPENDIX

**Listing of a Short MATLAB Program That Uses the EYELINK
and Psychophysics Toolboxes to Create a Simple Gaze-Dependent Display**

```
% Short MATLAB example program that uses the EYELINK and Psychophysics
% Toolboxes to create a real-time gaze-dependent display.

% STEP 1
% Initialization of the connection with the EYELINK Gazetracker.
% exit program if this fails.

if (EYELINK('initialize') ~= 0)
    return;
end;

% STEP 2
% Open a graphics window on the main screen
% using the PsychToolbox's SCREEN function.

screennr = 0; % use main screen
[window, screenRect]=SCREEN(screennr, 'OpenWindow', 0);
white=WhiteIndex(window);
black=BlackIndex(window);

% STEP 3
% Provide EYELINK with details about the graphics environment
% and perform some initializations. The information is returned
% in a structure that also contains useful defaults
% and control codes (e.g. tracker state bit and EYELINK key values).

el=initeyelin defaults;

% make sure that we get gaze data from the EYELINK

EYELINK('command', 'link_sample_data = LEFT,RIGHT,GAZE,AREA');

% open file to record data to

EYELINK('openfile', 'demo.edf');

% STEP 4
% Calibrate the eye tracker using the standard calibration routines
EYELINK('trackersetup');

% do a final check of calibration using driftcorrection
EYELINK('dodrftcorrect');

% STEP 5
% start recording eye position
EYELINK('startrecording');

% record a few samples before we actually start displaying
waitsecs(0.1);
% mark zero-plot time in data file
EYELINK('message', 'SYNCTIME');
stopkey=KbName('space');
eye_used = -1;
```

APPENDIX (Continued)

```

% STEP 6
% show gaze-dependent display
while 1 % loop till error or space bar is pressed
    % Check recording status, stop display if error
    error=EYELINK('checkrecording');
    if(error~=0)
        break;
    end
    % check for keyboard press
    [keyIsDown,secs,keyCode] = KbCheck;
    % if spacebar was pressed stop display
    if keyCode(stopkey)
        break;
    end
    % check for presence of a new sample update
    if EYELINK('newfloatsampleavailable') > 0
        % get the sample in the form of an event structure
        evt = EYELINK('newestfloatsample');
        if eye_used ~= -1 % do we know which eye to use yet?
            % if we do, get current gaze position from sample
            x = evt.gx(eye_used+1); % +1 as we're accessing MATLAB array
            y = evt.gy(eye_used+1);
            % do we have valid data and is the pupil visible?
            if x~=el.MISSING_DATA & y~=el.MISSING_DATA & evt.pa(eye_used+1)>0
                % if data is valid, draw a circle on the screen at current gaze position
                % using PsychToolbox's SCREEN function
                gazeRect=[ x-7 y-7 x+8 y+8];
                SCREEN(window, 'FrameOval', white,gazeRect,6,6);
            else
                % if data is invalid (e.g. during a blink), clear display
                SCREEN(window, 'FillRect',black);
            end
        else % if we don't, first find eye that's being tracked
            eye_used = EYELINK('eyeavailable'); % get eye that's tracked
            if eye_used == el.BINOCULAR; % if both eyes are tracked
                eye_used = el.LEFT_EYE; % use left eye
            end
        end
    end
    end % if sample available
end % main loop

% wait a while to record a few more samples
waitsecs(0.1);

% STEP 7
% finish up: stop recording eye-movements,
% close graphics window, close data file and shut down tracker

EYELINK('stoprecording');
SCREEN(window,'close');
EYELINK('closefile');
EYELINK('shutdown');

```

(Manuscript received March 20, 2002;
revision accepted for publication July 27, 2002.)