

INDEX

A

- Affinity tag
cellulose-binding module (CBM).....35–44
chitin binding domain (CBD).....48
cleavage
 autoprocessing59–69
 intein..... 41, 44, 48, 49, 53, 54, 56, 57
 tobacco etch virus (TEV) 115, 123
1D4 epitope tag 2, 3, 5–7, 10
elution..... 14, 22, 27, 30, 48, 78, 108, 115, 122, 145, 198
flag tag..... 2, 4, 138, 145, 148, 151, 196
glutathione-*S*-transferase (GST) tag 59, 72, 121,
 195, 196, 198–202
green fluorescent protein (GFP) tag..... 38–41, 43,
 121, 151–160, 164, 172, 176, 179
hemagglutinin (HA)..... 2, 60–61, 121, 136, 196, 200
heme-tag.....17–32
hexahistidine-tag (His-tag) 2, 59, 107, 108, 115, 118
his-MBP tag (Maltose Binding Protein).....72,
 81–93, 116
maltose binding protein (MBP) 72, 81–93, 116
multifunctional GFP (mfGFP) 152, 156, 158
Myc 2, 152
small ubiquitin-related modifier (SUMO) 71–79, 194
strep-tag.....107–108, 116, 117, 158, 196
streptavidin-binding peptide (SBP)..... 122, 125–127,
 152, 156–158, 176, 200
synaptobrevin Avi purification-tag (SnAvi).....163–173
tamavidin.....95–105
- Agarose
 beads.....64, 66, 69, 122, 124, 126, 130, 131,
 136, 139, 141, 143–145, 148, 201, 203, 204
 gel50, 51, 63, 64, 82, 86, 99, 179, 180, 182
- Antibiotic resistance marker
 ampicillin..... 21, 48, 50, 51, 63, 65, 67,
 83–87, 93, 100, 110, 112
 blastidicin 139, 141, 153, 155
 chloramphenicol 20, 21, 84, 87, 92
 Hygromycin.....139, 140, 153, 155,
 159, 178, 179, 182, 183, 187, 189
 kanamycin.....77, 93
 Puromycin 178, 182, 187
- Antibody
 monoclonal.....2–7, 13, 179
 polyclonal.....1, 3, 13

B

- Beads
 activated.....8, 26
 agarose64, 66, 69, 122–124,
 126, 130–132, 136, 139, 141, 143–145, 148, 184,
 201, 203, 204
 anti-FLAG M2 139, 145, 184
 biotinylated102, 136, 137, 147, 172
 dyna- 98, 102
 glutathione..... 201, 203, 204
 IgG 122, 123, 126, 131, 132, 145
 magnetic8
 micro- 96, 98, 102–105
 nickel NTA205
 protein A168
 protein G 122, 123, 126, 131, 139, 143, 145, 148
 sepharose 8–10, 13, 18, 20, 26, 31
 Strep-Tactin..... 154, 156–158
 streptavidin 122–124, 126, 127, 131, 132,
 136, 139, 141, 144, 172, 184–186, 205
- Biotin
 biotinylated96, 98, 102, 136, 163, 171
 in vivo biotinylation..... 135–148, 171
- Blotting
 immuno- 101, 156, 157, 207
 western2–4, 7, 10–13, 108,
 114, 117, 124, 125, 127, 130–132, 145, 146, 148,
 166, 170, 172, 176, 179, 182, 183, 185, 197, 201,
 206, 207
- C
- Cells
 bacterial 1, 5, 60, 68, 72, 87, 100, 160, 171
 BL21..... 24, 63, 65, 68, 77, 110, 112, 170
 competent.....49, 51, 55, 63, 65, 82, 83, 86, 92,
 110, 177, 181
 DH5 α 51, 55, 63, 65, 83, 86, 92, 93, 177
 Escherichia coli.....71, 82, 83, 85–93, 97
 eukaryotic 60, 72, 122, 125, 156, 164
 gyrA⁺..... 83, 85, 92
 HEK293T3–6, 9–10, 12, 153,
 155–156, 166, 169, 198, 202
 HeLa166, 169
 mammalian 1, 5, 153, 155,
 164–166, 168, 170, 171, 175–190

Cells (cont.)

- mouse embryonic stem cells (mESC) 135–148
- sorting (FACS) 176, 179, 182, 183
- T1^R 82, 85
- T-REx 155, 159
- U2OS 182–186, 189
- yeast 1, 5, 36, 122, 175, 176

Chromatography

column

- chitin 50, 52, 55
- glutathione (GSH) 203
- HIS 27, 29, 118, 157
- Ni-NTA 114, 119, 154
- RAC 40, 42
- strep-tactin 154, 157, 158

- elution 25–26, 29, 35, 40, 42, 77, 78, 108, 115, 116, 118, 119, 157, 158

- immobilized metal affinity chromatography (IMAC) 60, 81, 107, 154, 157, 158

Cloning

- deoxyribonucleic acid (DNA) 49–51, 55, 61, 63–65, 75, 82, 83, 85, 86, 98–100, 154, 159, 177, 180, 181
- plasmid 20, 49, 50, 63, 65, 83, 85, 86, 125, 154–155, 178, 181
- polymerase chain reaction (PCR) 5, 48–50, 56, 64, 65, 75, 77, 82, 83, 85, 86, 93, 98, 154, 155, 180–182
- recombinant cloning 82–86, 96, 98–100, 154, 155, 159
- restriction enzyme 49–51, 64, 154, 155, 177
- transfection 182
- transformation 51, 55, 86, 93, 99, 181, 183
- vector 20, 49–51, 53, 60, 63–65, 67, 75, 77, 82–86, 93, 98, 99, 155, 177, 178, 180–182, 189, 202

D

Detection and confirmation

- co-immunoprecipitation 2, 3
- coomassie 3, 20, 42, 69, 97, 157
- gel electrophoresis 97, 157
- GFP fluorescence 152, 157, 158
- in vivo 71, 72, 74, 78, 138, 207
- LC-MS 177, 187, 189
- mass spectrometry 78, 125, 172
- proteomics 136
- SDS-PAGE 20, 42, 43, 69, 78, 97, 100–102, 117, 118, 129, 157, 158
- western blot 3, 7, 117, 125, 166, 170

E

- Epitope 1–14, 91, 121, 136, 152, 164, 166, 176

I

- Immunoaffinity 3, 5, 8–9, 11, 13, 14
- Inducible
 - expression 24, 59–69, 128, 153, 155–156
 - promotor 24, 48, 128
- Intein 36, 38, 41–44, 48, 49, 53, 54, 56, 57

L

- LC-MS/MS 176, 177, 184, 186–189, 201
- Lysis
 - buffer 25, 38, 41, 64, 66, 68, 76, 77, 123, 125–127, 130, 160, 165, 167, 168, 170, 179, 180, 184–186, 190, 196, 199, 200, 202, 203
 - cell lysis 38, 41, 52, 56, 68, 75, 76, 79, 165, 167, 168, 196

M

- Mass spectrometry 11, 78, 121, 124, 125, 127, 128, 132, 136, 148, 171, 172, 177, 184–186, 193, 196, 206–207
- Membrane
 - blot 2–4, 7, 10–13, 100, 102, 108, 117, 129, 166, 179, 183
 - cell 1–6, 9, 10, 12, 13, 110–111, 113, 183
 - extract 2, 9–10
 - protein 1–7, 9–14, 18, 26, 31, 97, 100–102, 110–111, 113, 116, 117, 119, 129, 166, 183
 - PVDF 7, 11, 97, 100, 102

P

- Peptide 3, 5, 6, 10, 11, 14, 17–20, 30, 36, 61, 73, 98, 103, 122, 125–127, 136, 139, 145, 148, 151, 152, 158, 163, 164, 170, 176, 179, 184–187, 196, 198, 206, 207
- Plasmid 3, 5, 9, 13, 20, 24, 38, 39, 41, 49, 50, 63, 65, 75, 83, 85, 86, 89, 92, 110, 112, 125, 128, 129, 138–140, 153–155, 166, 170, 178, 181, 187, 189, 200, 202
- Polymerase chain reaction (PCR) 5, 48–50, 55, 56, 63–65, 68, 75, 77, 82, 83, 85, 86, 91, 93, 98, 99, 154, 155, 180–182, 207
- Protease 2, 6, 9, 59–61, 64, 68, 69, 72, 73, 84–91, 111–116, 118, 119, 122–124, 126, 130, 131, 139, 143, 144, 147, 154, 157, 160, 164–166, 169, 170, 172, 176, 179, 183, 184, 186
- Protease inhibitor 6, 9, 64, 68, 111, 113, 114, 123, 124, 130, 139, 143, 144, 154, 160, 165, 169, 179, 183, 184
- Protein
 - bait 121–125, 127–129, 131, 132, 187, 193, 194, 202, 203, 207
 - complex 2, 36, 135–148, 152, 153, 156, 157, 163, 164, 169, 171, 175, 184, 186, 193–196, 204

elution.....10, 11, 14, 22, 27,
 43, 44, 48, 3, 54, 62, 64, 76, 79, 87, 90, 108, 118,
 122, 124, 127, 129, 132, 148, 154, 156, 158, 170,
 186, 199, 204

endogenous.....136, 138, 141, 147,
 170, 171, 190, 197, 198

expression 1, 3, 11, 13, 17, 19–21, 24, 25,
 30, 31, 36–41, 49–52, 61, 65, 67, 68, 72, 73, 75,
 77–79, 81–93, 95–105, 112, 113, 117, 138,
 151–153, 155, 156, 164, 170, 171, 176, 180, 182,
 183, 189, 190, 197, 200, 202, 207

extraction22, 51, 84, 87, 100, 177

folding68, 71, 72, 74, 91, 116

fusion..... 17, 48–52, 60, 61, 64, 67, 69,
 71–79, 81, 82, 84–93, 95–105, 110–112, 116,
 118, 121, 151–153, 155–158, 160, 164, 165,
 76, 183, 184, 186, 202

interaction.....2, 95, 107–108, 117,
 121–132, 151, 152, 175, 177, 188, 190, 193

localization.....2, 7, 12–13, 151, 152, 164, 176, 177, 183

membrane.....1–14, 18, 97, 100–102,
 110, 111, 113, 116, 117, 119, 129, 166, 183

modification177

partner 2, 72, 74, 96, 98–101, 103–105,
 108, 116, 121–124, 128, 137, 138, 169, 171, 176,
 185, 186, 188, 190, 193, 194, 196

production 24, 35, 61, 81, 85, 86, 91, 96, 117

proteomics136

purification121–132
 bimolecular affinity purification (BAP)193–208
 multifunctional affinity purification
 (MAP)175–190
 tandem affinity purification (TAP)..... 121–132,
 136, 152, 153, 157, 163, 164, 166–169, 175,
 193–208

quantification.....5, 13, 20, 23, 26,
 29–30, 61, 91, 122, 129, 130

recombinant.....35–44, 47, 59, 61, 62, 71,
 81, 85, 95–105, 107, 108, 112, 116, 117,
 151–160, 170, 195, 196, 204, 207

solubility24, 31, 59, 61, 62, 68, 69,
 71–79, 81, 84, 86–89, 91, 93

stability 35, 59, 61, 204

ubiquitinated..... 196–198, 202, 204, 205

R

Restriction

digestion 51, 64, 65, 67, 68, 99, 181, 206
 enzyme 38, 49–51, 64, 99, 103, 154–155, 177
 site 5, 61, 67, 98, 99, 154–155, 177, 180

S

SDS-PAGE..... 7, 10, 11, 20, 22–23, 28–29, 39–43,
 52, 53, 56, 57, 62, 63, 65–67, 69, 77–79, 84,
 87–89, 93, 97, 100–104, 108, 114, 115, 117, 118,
 124, 127–129, 131, 132, 140, 142, 146–148, 157,
 158, 166, 169, 172, 179, 183, 196, 205, 206

Solubility

enhancer2, 81, 89
 protein 1–14, 24, 31, 59, 71–79, 86,
 88–91, 95, 96, 98, 100–103, 105, 108, 118, 119,
 128, 130, 160

tag..... 1–14, 24, 59, 61, 62, 68, 69,
 71–79, 81, 84, 86–91, 93, 95, 96, 108, 129

Supernatant 8, 10, 25, 38–42, 52, 57, 60, 62,
 65–67, 69, 77–79, 88, 100, 102, 103, 113, 114,
 122, 125–127, 131, 141–145, 148, 156, 158, 160,
 167–172, 184, 203, 205

T

Tag

activation 6, 8, 9, 22, 26, 31, 43, 48, 53,
 59, 60, 68, 73, 108, 110, 117, 119, 163, 171, 196
 affinity2–5, 13, 14, 17–32, 35–44,
 47–57, 59, 60, 72, 74, 78, 81, 95, 96, 107, 108,
 112, 115–118, 121–132

epitope1–14, 91, 121, 136, 152, 164, 166, 176
 purification1–14, 17–32, 35–44, 47–57,
 59–69, 71, 72, 74, 76–79, 84, 86–91, 95, 96, 107,
 108, 112, 115–118, 121–131, 136–138, 140,
 151–160, 163–173, 175–178, 182–187, 190,
 194–202, 207

removal 8–12, 14, 18, 19, 25, 26, 28, 29,
 31, 36, 38–40, 42, 43, 47, 48, 59, 65–67, 69,
 71–74, 79, 88, 90, 93, 102, 108, 115, 118, 123,
 125–127, 131, 176, 186

solubility 1–14, 24, 59, 61, 62, 68, 69,
 71–79, 81, 84, 86–91, 93, 95, 96, 108, 129