

# Index

## A

- Adverse reaction to metal debris (ARMD), 10
- ALVAL. *See* Aseptic lymphocyte dominated vasculitis associated lesions (ALVAL)
- Articular surface replacement (ASR)
  - arc angles, 81
  - cup inclination, 81
  - lower blood metal, 81
  - XL systems, 74
- Articulation related failure
  - catastrophic wear, 10
  - destroyed ceramic head, 10
  - low wear bearing surfaces, 9
  - malposition, 15
  - materials, 9
  - osteolysis and head eccentricity, 10
  - revision
    - ceramic on ceramic, 13–14
    - decision-making factors, 11–12
    - metal on metal, 13
    - PE, 12
    - treatment, 11
  - X-PE liners, 10
  - x-rays, 14, 15
- Aseptic lymphocyte dominated vasculitis
  - associated lesions (ALVAL), 10–11
- ASR. *See* Articular surface replacement (ASR)

## B

- BIOLOX® delta, 127, 129
  - chromium, 146
  - crack extension, 147
  - description, 146
  - fracture toughness and strength, 147
  - microstructure, 146, 147
  - physical properties, 146
  - reinforcing mechanism, 147, 148
  - zirconia particles, 148
- BioloX Delta alumina matrix composite
  - bearings, 5
- BioloX Forte alumina ceramic-on-ceramic
  - bearings, 5
- Birmingham hip resurfacing (BHR), 81

## C

- Ceramic composite, long-term stability
  - BIOLOX® delta, 146–148
  - component and material testing
    - ball head fracture, 148, 149
    - burst tests, 148
    - stress levels, 149
  - degradation, 146
  - description, 145
  - lifetime experiments
    - autoclaving, 151, 152
    - bending tests, 149, 150
    - inherent protection, 152
    - monoclinic phase, cyclic loading, 151, 152
    - phase transformation, 151
    - residual strength and monoclinic phase, 151
    - test matrix, 150, 151
  - lifetime limiting effects
    - fatigue, 150
    - human body environment, 150
    - hydrothermal, 150
    - long-term treatment, 149
- Ceramic-on-ceramic bearings
  - BioloX delta, 127–128
  - clinical review
    - alumina, 157
    - aseptic loosening, 156
    - debris, 155
    - fixation, 155
    - HA bonding, 156
    - Harris hip scores, 156
    - methods, 156
    - plastic debris, 157
    - prosthetic implants, 155
  - delta ceramics, 141
  - demographic analysis, 138
  - description, 127
  - fracture rate
    - alumina, 130
    - femoral heads, 130
    - in vivo*, 130, 131
    - osteolysis, 131

- Ceramic-on-ceramic bearings (*cont.*)
- handlings
    - intraoperative damage, 140
    - shoehorn, 140
    - suction device, 139
  - hip noise, 138–139
  - implant failure, 127, 128
  - laser etching, 127
  - materials and method
    - prosthesis system, 136
    - retrospective analysis, 136
    - Zimmer™ Alloclassic Variall® system, 136, 137
  - noises, 137
  - pelvic protrusion, 137
  - properties and burst strength, 128, 129
  - reputation, 127, 128
  - revision
    - breakage, 13–14
    - stable cup, 13
    - strategy, 133
  - risk factors, 131–132
  - squeaking, 133–136
  - stripe wear, retrieval liner and head, 138
  - survivorship, 129–130
  - transient squeaking, 141
  - tribological properties, 130
  - wear rate, 5
- Ceramic-on-metal bearings, 4–5
- Conventional vs. sequentially cross-linked polyethylene liners
- bone resorption, 33
  - larger femoral head sizes, 33–34, 36
  - materials and methods
    - liners and femoral heads, 34, 35
    - six-station hip wear simulator, 34
    - wear testing, 34
  - results
    - average gravimetric wear, 35, 36
    - X3T liners wear, 35
  - simulator studies, 37
- Cross-linked polyethylene (X-PE)
- breakage, 10
  - failures, 36–37
  - larger femoral heads, 34, 36
  - THR, 33
  - wear and osteolysis reduction
    - acetabular osteolysis, 65, 66, 68
    - annual wear rates, 65
    - duration vs. conventional PE, 67–68
    - HXLPE, 68
    - patient characteristics and implant parameters, 64
    - patients and methods, 61–63
    - randomization, 64
    - RCTs, 66
    - total linear head penetration vs. duration and conventional PE, 64–65
    - wear rates, 37, 67
- D**
- Delayed-type hypersensitivity (DTH)
- CoCr debris, 113
  - Sikomet articulation and ceramic-on-poly-ethylene, 122
  - T-cells, 122
- Durom. *See* Large-diameter metal-on-metal cup (LDMOM) and standard stem
- E**
- ECiMa™
- annealing technique, 170
  - antioxidants, 177
  - average oxidation, samples, 175–177
  - fatigue crack propagation, 173–175
  - free radicals, 170
  - HXLPE, 178
  - ionising radiation, 169
  - mechanical properties, 173
  - methods
    - average oxidation, samples, 172–173
    - fatigue crack propagation, 171
    - mechanical properties, 170–171
    - oxidation, 172
    - wear testing, 171–172
  - oxidation, 175, 176
  - peri-prosthetic osteolysis, 169
  - UHMWPE, 177
  - vitamin E, 178, 179
  - wear testing, 174, 175
- Elongation at break (EAB), 24, 25
- F**
- Free radicals
- crystalline region, 169
  - HXLPE materials, 169, 170
  - radiation stabilisation, 170
- H**
- Harris hip scores (HHS), 62, 64, 156
- Highly cross-linked polyethylene (HXLPE)
- fatigue crack, 173, 174
  - ionising radiation, 169
  - mechanical properties, 173
  - stress intensity, 175
  - vitamin-doped, 27

- Hip prostheses  
 ceramic-on-ceramic bearings, 5  
 ceramic-on-metal bearings, 4–5  
 metal-on-metal bearings, 4  
 polyethylene bearings  
   head size, 3  
   osteolysis, 5  
   wear rate, 4
- HXLPE. *See* Highly cross-linked polyethylene (HXLPE)
- HXLPE isoelastic monoblock pressfit cup, oxidation prevention  
 ceramic-on-ceramic bearings, 22  
 clinical documentation, 25  
 hip simulator testing  
   inlays, 24  
   materials, 24  
   test fluid, 24  
 mechanical properties  
   accelerated ageing, 25, 26  
   tensile testing, 25, 26  
 mechanical testing, 23–24  
 RM Pressfit vitamys cup, 23  
 UHMWPE and HXLPE, 22, 23  
 vitamin-doped HXLPE material, 22  
 vitamys material, 22  
 wear rates  
   accelerated ageing, 27  
   clinical follow-up, 27–29  
   cross-linking, 26–27  
   head material and diameter, 27
- I**
- Immunological adverse reaction  
 antibodies, immunohistochemical evaluation, 115, 116  
 Bicon cups, 120  
 CD3 and CD68, 121–122  
 ceramic-on-polyethylene bearing, 116  
 CoCr debris, 113  
 DTH, 122  
 fibrin exudation, neo-capsule, 121  
 IV-delayed hypersensitivity, 120  
 metal hypersensitivity, 116  
 Milosev's series, 121  
 MoM and CP revision, 115, 117–118  
 MoM THAs, 113  
 necrosis, synovial membrane, 115, 119  
 patients  
   cementless Bicon-Plus, 114  
   ceramic femoral head, CP, 114  
   ChemMate™ Dako EnVision™, 115  
   CP and MoM, 114  
   giant cells and neutrophils, 115  
   lymphocytic infiltration, 115  
   periprosthetic tissues, 114  
   primary implantation, THA, 114  
   periprosthetic osteolysis, 113  
   perivascular infiltration, lymphocytes, 120  
   perivascular lymphocytic infiltration, 115, 119  
   Sikomet articulation, 121  
   Sikomet bearing, 119  
   T and B-lymphocytes, 115, 120  
   T-suppressor, 115, 121
- K**
- Kaplan-Meier (KM) method, 102
- L**
- Large-diameter metal-on-metal cup (LDMOM) and standard stem  
 cemented stem, 92  
 clinical assessment, 92  
 data collection, 93  
 Durom acetabular component, 91  
 LHMOM, 91  
 longevity and pseudotumors, 94  
 low threshold patients  
   acetabular cup, 96  
   cup loosening, 95  
   loose stems, 95  
   metallosis, 95  
 materials and methods, 91–92  
 opening angle, 91–92  
 osteoarthritis (OA), 94  
 prosthesis infection, 94  
 radiographic analysis, 94  
 radiological analysis  
   anteroposterior pelvic and lateral radiographs, 92  
   cup inclination, 92  
   heterotopic ossification, 93  
   loosening definition, 92  
   radiographs, 92–93  
 revision rate, 95  
 Trendelenburg sign, 93
- M**
- Metal-backed cup, 99–100  
 Metal-on-metal bearings  
   acetabular cup, 73  
   femoral bone stock, 73  
   imaging, soft tissues, 85  
   immunological adverse reaction (*see* Immunological adverse reaction)  
   measuring component position, 84–85

- Metal-on-metal bearings (*cont.*)  
 MHRA and ASR, 74  
 pathogenesis, 73–74  
 revision, 13  
 soft tissue reactions  
   pseudotumour, 86–87  
   revision surgery, 87–88  
 wear and ion levels  
   chromium and cobalt, 74, 75  
   component design, 81–82  
   diametric size tolerance, 76  
   high bearing surface, 77–78  
   linear/volumetric measurement, 76  
   material loss, large diameter, 83–84  
   parts per billion (ppb) and micrograms  
     per litre ( $\mu\text{g/L}$ ), 74  
   patient factors, 83  
   roundness measurement and CMM, 76  
   sensitivity and specificity, 76  
   soft tissue pathology, 74, 75  
   surgical positioning (*see* Surgical  
     positioning)  
   wear rate, 4  
 Milosev's series, 121  
 Mooring screw, 163, 165
- N**  
 Novae® acetabular cup, 165
- O**  
 Orthopaedic arthroplasty. *See* Vitamin E-  
 blended UHMWPE
- P**  
 Perkin-Elmer DSC6 differential scanning  
 calorimeter, 43  
 Perkin-Elmer SPECTRUM SPOTLIGHT 300  
 auto-image microscope, 44  
 Polyethylene (PE) bearings  
   head size, 3  
   revision, 12  
   wear rate, 4  
 Polyethylene (PE) wear  
   anteroposterior (AP) radiograph, 101, 102  
   aseptic loosening, 106  
   cemented hip prostheses, 99  
   conventional UHMWPE, 107–108  
   femoral head, 101  
   head penetration rate, 101  
   heterogeneity, 106  
   highly cross-linked polyethylene  
     (HXLPE), 105  
   hydroxyapatite (HA), 99  
   implant failure  
     KM survival estimation, 104, 105  
     liner and related revision, 103, 104  
     index surgery and UHMWPE liners, 105  
     Kaplan-Meier (KM) method, 102  
     macrophages, 99, 100  
     metal-backed cup, 99, 100  
     metal-on-metal (MoM) bearing, 106  
     osteolysis  
       boxplot, 102, 103  
       cup inclination angle, 103  
       mean measurement, 103, 104  
     patient demographics, 101  
     production and implantation, 100  
     radiograph, 106  
     squeaking, 106  
     steeper inclination angle, 106  
     total hip prostheses, 101  
     treatment algorithm, 106, 107  
     uncemented hip prostheses, 100  
 Postel Merle d'Aubigné score, 163  
 Pseudotumour, 86–87
- R**  
 Reicher-Jung POLY-CUT S microtome, 44  
 RM Pressfit vitamys implantation, 23, 25  
 Roman V 1.70 software, linear head  
 penetration, 63
- S**  
 Sikomet articulation, 121  
 Sikomet bearing, 119  
 Soft tissue reactions  
   pseudotumour, 86–87  
   revision surgery, 87–88  
 Squeaking  
   ceramic fractures, 133  
   incidence, 133, 134  
   influencing factors  
     component orientation, 135  
     fluid film disruption, 135  
     patient-related factors, 135–136  
     prosthetic design and material, 134–135  
 StatView statistical software, 163  
 Surgical positioning  
   articular arc angle, 80  
   cup inclination and negative correlation, 78, 79  
   femoral diameter, 79  
   radius, contact patch, 80  
   wear and metal ion, 78
- T**  
 Titanium dual-mobility socket  
   cementless fixation, 166  
   femoral stem loosening or migration, 164

- intra-prosthetic dislocation, 164, 166  
 loosening, acetabular cup, 163  
 low friction coefficient, 167–168  
 material and method, patients  
   clinical and radiographic  
   examination, 163  
   etiologies, 162–163  
   mobile liner, 163  
   primary hip arthroplasty, 161  
   StatView statistical software, 163  
   titanium Novae, 162  
 mooring screw, 163, 165  
 Novae® acetabular cup, 165  
 osteolysis, 163  
 polyethylene component, 166  
 Postel Merle d' Aubigné score, 163  
 prosthetic instability, 161, 167  
 radiography, THA, 163, 164  
 wear, 167  
 Trendelenburg sign, 93
- U**
- UHMWPE. *See* Ultra-high molecular weight  
 polyethylene (UHMWPE)  
 Ultra-high molecular weight polyethylene  
 (UHMWPE). *See also* Vitamin  
 E-blended UHMWPE  
 acetabular components, 169  
 cross-linking, 33  
 fatigue strength, 177  
 implants  
   accelerated ageing, 24  
   gamma-sterilised, 26  
   lifetime, 23  
   mechanical properties, 26  
   wear rate, 27, 28  
 plastic deformation, 169  
 wear resistance, 60  
 Uncemented prostheses. *See* Polyethylene  
 (PE) wear
- V**
- Visual analogue scale (VAS), 28  
 Vitamin E-blended UHMWPE  
   biocompatibility, *in vitro* and *in vivo*  
   studies  
   machining, 44  
   results, 50–51  
   tests and conditions, 44, 45  
 chemico-physico-mechanical  
 characterisation  
   concentration profile, 49  
   crystallinity, samples, 43, 46
- FTIR measurements, 44, 46  
 IR spectra, accelerated ageing test, 48  
 non-additivated material, 46, 47  
 oxidation levels, 47, 48  
 quantitative amount, ketones  
   concentration, 48, 49  
 Reicher-Jung POLY-CUT S  
   microtome, 44  
   stress-strain curve, 50  
   tensile test data, 50  
 cumulative mass loss, 55  
 DSC measurements, 53  
 g-irradiation, 41  
 GUR, 42  
*in vivo* mechanical stress and oxidation, 42  
 post-irradiation treatment elimination, 54  
 samples, disposal, 43  
 wear test  
   behaviour determination, 44  
   contact areas, 52  
   EDTA, 45  
   liner areas, microscope, 45, 46  
   rate, non-additivated UHMWPE, 51  
   servo hydraulic six-station hip  
   simulator, 45  
   surface damage, 52–54
- W**
- Wear and osteolysis reduction, X-PE  
   acetabular osteolysis, 68  
   duration vs. conventional PE, 67–68  
   HXLPE, 68  
   patients and methods  
     ABG-II implant system, 61  
     assessment, 62  
     conventional and duration PE, 62  
     Fisher's exact test, 63  
     linear head penetration measurement, 63  
     polyethylene inserts, 62  
     unilateral THA, 61  
 RCTs, 66  
 results  
   acetabular osteolysis, 65, 66  
   annual wear rates, 65  
   patient characteristics and implant  
   parameters, 64  
   randomization, 64  
   total linear head penetration vs. duration  
   and conventional PE, 64–65  
   wear rates, 67
- Z**
- Zimmer™ Alloclassic Variall® system, 136, 137